

Non-Invasive Miniature Bioreactor for Laboratory Engineering **(NIMBLE)**


Project Sponsor:

Todd Thorson and Lincoln Labs

Mentors:

Douglas Hart, Juergen Schoenstein, Jared Berezin,
Rohan Sanghai, Jessica Lam

05 December 2024



Miguel Talamantez • Timber Carey • Henry Smith • Christian Duessel • Rakibul Chowdhury
Sehar Lohawala • Reyna Ayala • Gerardo Berlanga • Faris Elnager • Abraham Lemaitre

What is a bioreactor?

- Vessel or tank that creates optimal environment for biological reactions
- Used for many applications
 - Cell culture
 - Fermentation
 - Biopharmaceutical production
- Come in many different sizes.



100mL-40L single-use bioreactors (Eppendorf)




250mL single-use bioreactor (Sartorius)



5L-30L reusable bioreactor (Sartorius)

Types of Bioreactors

- Batch
 - All nutrients added at the start of the process
- Fed-batch
 - Nutrients are added
 - No products are removed
- Continuous
 - Products are simultaneously removed
 - Higher population concentrations
 - Remove waste and product during reaction
 - Continuously add nutrients



x = culture population v = culture volume s = substrate concentration


How NIMBLE works

Set Up

- Add media to vessel
- Assemble bioreactor
- Calibrate sensors
- Attach tubing
- Attach feed, acid and base inputs



Input control parameters:

- Feed rates
- Target temp.
- Target pH
- Stir rate




How NIMBLE works

Run Reaction



How NIMBLE works

Clean Up



Our Niche




	Our Bioreactor	General Large-Scale	General Lab-Scale
Low-cost (reuse)	✓	□	□
Single Use	✓	□	□
Temp sensing and control	✓	✓	✓
pH sensing and control	✓	✓	□
DO sensing and control	✓	✓	□
OD sensing	✓	✓	✓
Perfusion Continuous Flow	✓	□	□
Chemostat Continuous Flow	✓	□	□
100% non-invasive sensing	✓	□	□
Lab-scale (~100mL)	✓	□	✓
Data Logging	✓	□	□

Have

Some have

Don't have

NIMBLE Assembly




Assembled Bioreactor

Disposable


Reusable

Reactor Vessel +
Stirring Motor


NIMBLE Assembly



8 Basic Design Requirements




5 User Requirements



1 Cost Requirement


Cost (Hard Req.)		
Requirement	Metric	Source
Overall cost	<\$1000	Team

Disposable Components




Innovative Vessel Body

- 100mL volume container
- Material: Nylon-12 (mass production)
 - Prototype: Formlabs biomed resin
- Trenches allow for non invasive sensing
- Baffled to ensure homogeneity of solution.





Vessel Lid

- Inputs and outputs for functional requirements
- Output draws from mid-vessel
- Internal mixer
 - Keyed interface for motor mounting
 - Motor coupled to stir bar through the lid
- 2 options (decided by user)
 - Option 1: PVDF membrane (filter) in the output
 - Option 2: no filter




Innovative Agitation Design




Example run at 300 RPM


Total turbulent kinetic energy before vs after baffles included



Oxygen Transport (at 300 RPM) without Gas Input




Reusable Components




Assembled Bioreactor

Reusable Components




Bioreactor Base




Reactor Vessel + Motor


Reusable Components




Stirring motor (not pictured)




Reusable Components



Stirring motor (not pictured)



Innovative 4-Sensor Setup



Sensing Platform Dimensions:
60x80mm


Dissolved Oxygen
(DO) Sensor

pH Sensing Channel

Optical Density (OD)
Sensor


Thermistor location

Sensor Package Exterior:
Bottom View




Reaction Chamber:
Side View

Dissolved Oxygen (DO) Sensor



Dissolved Oxygen
(DO) Sensor




Sensor Package Exterior:
Bottom View

Reaction Chamber:
Side View

DO Sensor


spec: $2-10 \text{ mg/L} \pm 0.2 \text{ mg/L}$

Dissolved Oxygen Sensing Schematic [2]




Similar to as demonstrated by [1]

DO Sensor Package




DO Sensor: Cross Section

DO Sensor Model




~1.44mW of light incident upon the photodiode




$0.35 \times 1.44 = 0.504\text{mA}$ of generated current in the photodiode

DO Amplifier Circuit


1. Converts current signal to photodiode to voltage signal
 2. Amplifies signal so easily detectable with microcontroller analog input



pH Sensing Channel



pH Sensing Channel





Reaction Chamber Exterior:
Bottom View

Reaction Chamber:
Side View

pH Sensor


spec: $4\text{-}9 \pm 0.1$



Absorption of an 8 μm thick PANI film [3]

Fundamental Mechanism:
pH Dependent Absorbance of Polyaniline

pH Sensor Design



Cross section


Sensing "Trench"

PANI Sensing Layer


Red LED

Sensing
photodiode

Reference
photodiode




Vessel Side View




pH Sensing Diagram


Optical Density (OD) Sensor



Reaction Chamber Exterior:
Bottom View





Reaction Chamber:
Side View



Typical Bacterial Growth
Curve


OD Sensor

spec: $0-2 \pm 0.1$ au



Similar to Pioreactor (competitor bioreactor)


OD Sensor Design



OD Sensor: Cross Section


Cost breakdown + Component Selection

Part	P/N	Supplier	Cost
LED	521-9704F	Mouser	\$0.67
Reference PD	SFH 213 FA	Mouser	\$0.60
2x Measuring PD	1621132	Jameco	\$0.95
2x amplifier circuits	Slide 19/20	Slide 19/20	Slide 19/20




OD Placement:
Bottom View

External Thermistor



Reaction Chamber Exterior:
Bottom View


Thermistor location




Reaction Chamber:
Side View

Thermistor Temperature Validity

spec: $25-50^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$




Thermal Network



Estimated sensing time constant: ~30 seconds

Temperature Sensing Requirement:
<30 seconds

Temperature Sensor Design



Thermistor Placement:
Bottom View

Cost breakdown +
Component Selection

Part	P/N	Supplier	Cost
Thermistor	NRCE154K 4200B1KS	Mouser	0.17
Thermal pad	B0BD5PJ6 FL	Suspirit	0.02

Heater Selection


12.7 x 2.54 cm 0.78W/cm² Flexible Silicone Rubber heat pad

Requirement	Metric
Temperature control	25°C-50°C ± 1°C
Optical density sensing	0-2 ± 0.1au
pH sensing	4-9
Leak and ingress proof	IP68

$$q = mc_w \Delta T \rightarrow 12,552J \approx 12kJ$$

- To heat to max in 10 min need **20W**
- Heat pad can supply 25W total
- Clamped near the bottom of the vessel

$$\frac{\text{Heated SA}}{\text{Volume}} \rightarrow \frac{32\text{cm}^2}{100\text{cm}^3} = 0.32\text{cm}^{-1}$$



Heater Thermal Insulation

Flexible Teknofibra 2mm Adhesive Sheet


Hard/Soft	Requirement
Hard	Ability to input parameters
Hard	Safety
Soft	Real-time display of parameters
Soft	Data logging
Soft	Portable

From Teknofibra Test Report:

- 300s Temp. Stabilization for $\Delta T = 275^\circ\text{C}$
- Assuming within 5% of final: $\tau = 100\text{s}$

$$T_{out} = T_{in} - \dot{Q}R_{wall}$$

Temperature @ External Wall: 38°C





Pump Selection

Kamoer KPMP10 Peristaltic Pump

Requirement	Metric
Temperature control	25°C-50°C ± 1°C
Size - Benchtop scale	100 mL
Continuous Feed (Soft Req.)	0.1 - 5 mL/min
Run time	2 weeks


- 1 Inflow, 1 Outflow, 2 pH Control (Acid/Base Input)
- Pump Flow Rates: 0.1-5.9 mL/min
- MED55 PVC Tubing 1.6 x 3.2 mm ID x OD




Software and Electronics

Controls Overview

- Required environment regulation
 - Temperature
 - pH
 - Dissolved Oxygen
 - Flow rate




Requirement	Metric
Temperature control	25°C-50°C ± 1°C
Optical density sensing	0-2 ± 0.1au
Dissolved oxygen sensing	2-10 mg/L
pH sensing	4-9
Continuous Feed (Soft Req.)	0.1 - 5 mL/min



Temperature Control


- Thermal model to understand dynamics
 - Governing equation from 1st law
 - Neglecting air convection
 - Vessel contacts heating pad

$$m_{water} * C_{p,water} * \frac{dT}{dt} - \frac{1}{R_{eff}} * T(t) = -\dot{W}$$




Temperature Control

- Open loop transfer function
 - First-order response
- Compensator design
 - Proportional-Integral control
 - Achieve heating time and no steady state error
- Input power step response graph
 - Time constant = 20 min
- Temp control scheme
 - High initial power
 - Decrease power until equilibrium



pH Control


- Closed-loop feedback on-off controller with a deadband
- Acid/base addition
 - Acid: HCl, pH 0-1
 - Base: NaOH, pH 13-14
- User inputs
 - Target pH
 - Tolerance
 - Cycle time
- Small fixed amount of acid/base added by commanding pump on/off quickly



[3]


Dissolved Oxygen Control

- Closed-loop feedback on-off controller
 - Constants vary from vessel to vessel
 - Characterization testing required
- Gas addition
 - Oxygen O₂
 - Nitrogen N₂
 - Carbon Dioxide CO₂
- User inputs
 - Dissolved oxygen quantity
- Small fixed amount of gas added
 - Measure then add more if needed



[6]

Control Circuitry



Circuit Function Summary

120V AC input from wall

5V DC Regulated output

12V DC Regulated output

PWM Actuator control

Analog optical sensing with signal filtering

Power Electronics

Requirements:

120V AC → 12V DC

12V DC → 5V DC

PWM Control (25W peak)

Analog Signal Input


Hardware Selection:

12V Wall Adapter

LM7805 Voltage Regulator

L298N Gate Driver


ADS1115 Multiplexer



AC/DC Converter

P/N: GST60A12-P1J

12V Component	Steady Power Consumption
Pumps	0.5 W
Heaters	11 W
Stirring	1 W



Power Electronics

Requirements:

120V AC → 12V DC

12V DC → 5V DC

PWM Control (25W peak)

Analog Signal Input


Hardware Selection:

12V Wall Adapter

LM7805 Voltage Regulator

L298N Gate Driver

ADS1115 Multiplexer




5V Voltage Regulator Circuit
P/N: LM7805

Components Utilizing 5V DC Source

Raspberry Pi VCC

Analog Input Multiplexer

Gate Driver Logic



Power Electronics


Requirements:

120V AC → 12V DC
12V DC → 5V DC

PWM Control (25W peak)
Analog Signal Input

Hardware Selection:

12V Wall Adapter
LM7805 Voltage Regulator
L298N Gate Driver
ADS1115 Multiplexer




H-Bridge Working Principle Diagram
P/N: L298N

Components Utilizing H-Bridge PWM Control

Pump Motors

Heater

Stirring Motor



Analog Input

Requirements:

120V AC → 12V DC

12V DC → 5V DC

PWM Control (25W peak)

Analog Signal Input

Hardware Selection:

12V Wall Adapter

LM7805 Voltage Regulator

L298N Gate Driver

MCP3008 Multiplexer


8x Sensors

DO photodiode


7x Amplifier Circuits

Amplifier Circuit
(Slide 19)

Temperature Sensor




Data Bus



Project Implementation and Risk Management

Design Failure Mode and Effects Analysis (DFMEA)



Lifecycle



Disposable

Run time requirement: 2 weeks/300 hours

Failure modes:

- Material degradation in acidic/basic conditions
- Stirring assembly wear

Reusable

Limiting Component: Motor lifetime

Pump motor lifetime: **6,000 hours [4]**

Stirring motor lifetime: 10,000 hours


Usage assumptions:

- Pumps running continuously
- Room temperature usage
- No extreme forces are being applied
- Assume perfect models for electrical components

Overall Budget Breakdown

Category	Cost
Electronics and Sensors	\$143.17
Heating and Fluids	\$557.01
Mechanical and Structural Components	\$57.00

Total Bioreactor Production Cost: \$757.18



Approximate Vessel (disposable) Cost

\$58.00


Next Steps and Potential Innovation

Next Steps:



- Develop solution to foaming
- Gas box add-on

Potential Innovation:


- 4 in 1 non-invasive sensor package
 - Temperature
 - Dissolved Oxygen
 - Optical Density
 - pH



System Recap



Disposable



Assembled Bioreactor

Bioreactor Base

Reactor Vessel +
Stirring Motor

Appendix

BOM: Sensing & Electronics

#	Category	Component Name	Quantity	Cost per	Total Cost
1	pH Sensing	LED	1		\$0.00
2	pH Sensing	Photodiode	2		\$0.00
3	OD Sensing	LED	1	\$0.67	\$0.67
4	OD Sensing	Photodiode	2	\$0.95	\$1.90
5	OD Sensing	Reference Photodiode	1	\$0.60	\$0.60
6	DO Sensing	LED	1		\$0.00
7	DO Sensing	Filter	1		\$0.00
8	DO Sensing	Photodiode	2		\$0.00
9	Power Electronics	Wall Power Adapter			\$0.00
10	Temp Sensing	Thermistor			\$0.00
11	Logic components	Raspberry Pi 5	1	\$60.00	\$60.00
12	Logic components	Raspberry Pi Touch Screen	1	\$60.00	\$60.00
13	Logic components	Sensor Routing PCB	1	\$20.00	\$20.00
	Total Cost:			\$143.17	

BOM: Heating & Fluids

#	Category	Component Name	Quantity	Cost per	Total Cost
1	Agitation	Impeller	1	\$1.00	\$1.00
2	Heating	Heating Pad	1	\$71.11	\$71.11
3	Heating	TIM (1-2mm thick, one side sticky)	0.009	\$104.00	\$0.90
4	Motors/Pumps	Feed Inlet	1	\$121.00	\$121.00
5	Motors/Pumps	Main Output	1	\$121.00	\$121.00
6	Motors/Pumps	Acid Input	1	\$121.00	\$121.00
7	Motors/Pumps	Base Input	1	\$121.00	\$121.00
8	Motors/Pumps	Stirring motor	1		\$0.00
		Total Cost:		\$557.01	

BOM: Structures

#	Category	Component Name	Quantity	Cost per	Total Cost
1	Vessel	Membrane	1	\$22.00	\$22.00
2	Bioreactor Structure	Check Valve	1	\$10.00	\$10.00
3	Bioreactor Structure	Tubing	1	\$15.00	\$15.00
4	Bioreactor Structure	Bolts n shit	1	\$10.00	\$10.00
Total Cost:			\$57.00		

DO BOM

Electronics ->

Ru-Sensing layer


Part	Amount	Supplier	Cost
TEOS Precursor	10's mL	Sigma-Aldrich	\$70/L
Ethanol	10's mL	TBD	small
DI water	10's mL	TBD	small
ruthenium(II) dichloride complex	1's mg	Medchem Express	\$118/gram

Part	Specification	Supplier	Cost
LED	620 nm	Thorlabs	\$7
Photodiode	350-1100 nm	Thorlabs	\$17
Edge Filter	640 nm Longpass filter	Hoya	\$32
Reference PD	SFH 213 FA	Mouser	\$0.60
2x Resistor 1	500 ohms	TBD	small
2x Capacitor 1	1nF	TBD	small
2x Opamp	OP07	TI	\$1.20

Control Transfer Function

$$H(s) = \frac{T(s)}{W(s)} = \frac{1}{m_{water} * C_{p,water} * s - \frac{1}{R_{eff}}}$$

$$m_{water} * C_{p,water} * \frac{dT}{dt} - \frac{1}{R_{eff}} * T(t) = -\dot{W}$$



Pump Lifetime




Technical Parameters

KPMP10-ST

Motor Type	Stepper motor
Flow Rate Range	0.1-5.9ml/min
Power	10 W
Control Mode	Stepper motor driver
Weight	About 145 g
Rated Voltage	24V
Tube Lifetime	Depending on actual working conditions and media
Motor Lifetime	6000 h

Sensor Calibration Procedure

1. Remove new vessel from packaging
2. Fill new vessel with 10 mL of known pH
3. pH calibrates offset current
4. Rinse well and add 10mL of DI water
5. Leave lid partially covered for 10 minutes
6. Calculate equilibrium_oxygen(temp, altitude)
7. DO calibrates normalized intensity to equilibrium_oxygen
8. Pour out DI water
9. Fill vessel with starting contents for reaction
10. OD calibrates starting value



References

[1]

<https://sci-hub.ru/https://www.sciencedirect.com/science/article/abs/pii/S0263224114004527>

[2] McEvoy, Aisling K. et al., Dissolved Oxygen Sensor Based on Fluorescence Quenching of Oxygen-sensitive Ruthenium Complexes Immobilized in Sol-Gel derived Porous Silica Coatings

[3] Mackenzie C. R. Denton, Natasha P. Murphy, Brenna Norton-Baker, Mauro Lua, Harrison Steel, and Gregg T. Beckham. Integration of pH Control into Chi.Bio Reactors and Demonstration with Small-Scale Enzymatic Poly(ethylene terephthalate) Hydrolysis. *Biochemistry* 2024, 63, 1599–1607.

[4] <https://www.kamoer.com/us/product/params.html?id=10140>

[5] Bolmanis, Emils, Konstantins Dubencovs, Arturs Suleiko, and Juris Vanags. 2023. "Model Predictive Control—A Stand Out among Competitors for Fed-Batch Fermentation Improvement" *Fermentation* 9, no. 3: 206.
<https://doi.org/10.3390/fermentation9030206>

[6] <https://control.com/textbook/closed-loop-control/onoff-control/>